













Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

DISPOSABLE EP ELECTRODES			
720 00-S-C	Neuroline Electrodes 2,000 Case (80 x 25-pack) - 1 Case		
720 00-S-B	Neuroline Electrodes 1,000 Box (40 x 25-pack)		
720 00-S-A	Neuroline Electrodes 500 Box (20 x 25-pack)		
720 00-S	Neuroline Electrodes 25/pack		
RA069A	NeuroPlus Electrode (Size: Small, 7/8" square) - 600 Box		
RA069	NeuroPlus Electrode (Size: Small, 7/8" square) - 60/pack		
RA070A	NeuroPlus Electrode (Size: Large, 7/8 x 1 3/8") - 600 Box		
RA070	NeuroPlus Electrode (Size: Large, 7/8 x 1 3/8") - 60/pack		
WA069	VivoTab Electrodes - 25 patient supply (25 packs of 4)		
WA070	VivoTab Electrodes - 100 patient supply (100 packs of 4)		
EP EAR TIPS (for ABR & ECochG)			
ER3-14A	Disposable Foam Ear Tips for ER Insert Earphone (regular, 13 mm - yellow) - 50/pack		
ER3-14B	Disposable Foam Ear Tips for ER Insert Earphone (baby, 10 mm - beige) - 50/pack		
ER3-14C	Disposable Foam Ear Tips for ER Insert Earphone (jumbo, 18 mm - yellow) - 24/pack		
ER3-14D	Infant Ear Tips for ER Insert Earphone (3.5 mm - red) - 20/pack		
ER3-14E	Infant Ear Tips for ER Insert Earphone (4 mm - red) - 20/pack		
TIPtrode® - ELECTRODE EAR TIPS for non-invasive ECochG – GOLD-FOILED			
ER3-26A	TIPtrode® Electrode Ear Tips for ECochG - to be used with ER3-60 (regular, 13 mm) - 20/pack		
ER3-26B	TIPtrode® Electrode Ear Tips for ECochG - to be used with ER3-60 (baby, 10 mm) - 20/pack		
OTHER DISPOSABLES and CONSUMABLES			
RA061	Nuprep Abrasive Skin Prepping Gel - 1 Box of 6 tubes (non-conductive, non-sterile, requires cleaning of residual on the skin with Alcohol Prep Pads)		
S360-1	Spectra 360 Electrode Gel - 2oz - Single (for Non-invasive ECochG Tiptrode® electrodes)		
RA238	Dynarex Electrode Skin Prep Pad 100/box		
















Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

OAE EAR TIPS			
KR004	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 4 mm, Colour: Red) - 100/pack		
KR005	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 5 mm, Colour: Blue) - 100/pack		
KR006	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 6 mm, Colour: Green) - 100/pack		
KR007	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 7 mm, Colour: Yellow) - 100/pack		
KR008	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 8 mm, Colour: Red) - 100/pack		
KR009	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 9 mm, Colour: Blue) - 100/pack		
KR010	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 10 mm, Colour: Green) - 100/pack		
KR011	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 11 mm, Colour: Yellow) - 100/pack		
KR012	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 12 mm, Colour: Red) - 100/pack		
KR013	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 13 mm, Colour: Blue) - 100/pack		
KR014	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 14 mm, Colour: Green) - 100/pack		
KR016	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 16 mm, Colour: Yellow) - 100/pack		
KR036	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 3-6 mm, Colour: Yellow) - 100/pack		
MF003	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 3-8mm, Colour: Blue) - 100/pack		
MF007	Single-use Ear Tips for OAE Probe / Impedance Tip Adapter (Size 7-11mm, Colour: Green) - 100/pack		
GRS-TS210	Vivosonic GP OAE Ear Tip Set in Clear Box 156 assorted single-use ear tips, 3-13 mm		
RA063	Screening Ear Tip Set in Clear Box - 60 assorted single-use ear tips		
A90-08-2-0	R4.8S Otodynamic DPOAE Probe Tip - 100/pack		
A90-16-2-0	R5.8B Otodynamic DPOAE Probe Tip - 100/pack		
A90-17-2-0	R6.8B Otodynamic DPOAE Probe Tip - 100/pack		
A90-25-2-0	R7M Otodynamic DPOAE Probe Tip - 100/pack		














Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

A90-09-2-0	R8M Otodynamic DPOAE Probe Tip - 100/pack		
A90-10-2-0	R9M Otodynamic DPOAE Probe Tip - 100/pack		
A90-11-2-0	R11M Otodynamic DPOAE Probe Tip - 100/pack		
A90-19-2-0	R13M Otodynamic DPOAE Probe Tip - 100/pack		
UG-TS	Otodynamics UG OAE Ear Tip Set in Clear Box		
Generation 2 (G2) SPARE PARTS			
FA014	Vivocheck - Device for Integrity™ system self-check		
WA045	VivoAmp		A90
WA059	VivoAmp Lead Cable (Black, Blue, Green, and Red)		
WA065	Alligator Lead Cable for VivoAmp (Black, Blue, Green, and Red)		
WA075	G2 CV-Amp with Gold Cup Lead Cables - EEG prefilter and amplifier with connector and calibration in EEPROM		A91
WA074	Gold Cup Lead Cable (Black, Blue, Green, and Red)		
WA032	G2 Single Channel Amplitrode® Patented In-situ EEG Amplifier with cable, connector, and calibration in EEPROM (DISCONTINUED -- Quantities are limited)		A81






Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

WA031	G2 Dual Channel Amplitrode® Patented In-situ EEG Amplifier with cable, connector, and calibration in EEPROM (DISCONTINUED -- Quantities are limited)		A82
WA055	G2 ER-3C Insert Earphone Set with cable, connector, and calibration in EEPROM		ER-3C-800
WA072	G2 ER2 Insert Earphone Set for 500 Hz to 8kHz ABR (stimulus level limits) with cable, connector and calibration in EEPROM		ER2
ER3-28K	ECochG Hardware Kit (includes: WA 110, ER3-26A, ER3-26B, S360-1)		
ER3-04	White Tube Nibbons for ER Insert Earphones - 10/pack		
ER3-05	Velcro Clips for ER Insert Earphones - 6/pack		
ER3-06(ABR)	Impedance Tip Adapter - use economy Eartips with ER Insert Earphones - Set of 1 Right (red) and 1 Left (blue) tube		
ER3-21(ABR)	Replacement Front Tubes for ER Insert Earphone - Set of 2 Right (red) and 2 Left (blue) tubes		
WA056	G2 TDH39 Audiometric Headphone		H-801
WA110	ECochG Lead Cable for VivoAmp (replaces ER3-60 Lead Cable)		
WA101	G2 Bone Conductor 10 Ohm with Stainless-Steel Headband, 72" Cord with straight plug (requires Integrity software version 8.11.3 and above)		B81
WA030	G2 Bone Conductor 10 Ohm with Stainless-Steel Headband, 72" Cord with straight plug		B71W
P-3333	B-71 Metal Headband for B-71 Bone Conductor		
RA148	B-71 Elastic Headband for B-71 Bone Conductor		P-1188
WA034	G2 Vivosonic General-use OAE Probe with cable, connector, and calibration in EEPROM		P81-GP

Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

100399	Vivosonic OAE Probe Holder		
HAL2573	OAE Probe Cleaning Kit - mini-brush, mini-pick, and mini-loop		
WA035	G2 Otodynamic General-use OAE Probe with cable, connector, 5 couplers and calibration in EEPROM		P81-UG
P40-UG-C	Otodynamic Test Cavity		
WA113	G2 Vivosonic GX OAE Probe with cable, connector, and calibration in EEPROM (for Integrity software version 8.11.3 and above)		P81-GX
RA252	G2 Vivosonic GX OAE Probe Nose (pack of 5)		
RA253	G2 Vivosonic GX OAE Probe Cavity		
RA045	DPC Probe Couplers (x10)		
RA108	Body & Lid for UGD probe		
VIVO-BATT-KIT	Battery Pack Kit (A/C Adapter, Charging Cradle, 2 Battery Packs)		
VIVO-BATT-PACKS	Battery Packs x 2 (colour may vary)		
100351	Battery Pack A/C Wall Adapter		
100410	Battery Pack Cradle Assembly		
RA254	Battery Charger for Individual Ni-MH Cells		

Integrity™ V500 Generation 2 (G2)
PARTS AND ACCESSORIES CATALOGUE
 (For XP and G1 Systems, see endnote *)

WA115	Individual Ni-MH Battery Cells, pack of 8		
V51-GS	Adjustable General-use Strap with trigger snaps for secure placing VivoLink™ on older children and adults		
V51-PS	Set of 2 Adjustable Pediatric Shoulder Straps with trigger snaps - for secure placement of VivoLink™ on toddlers and young children		
100018	Printer Labels - 2 rolls of 300 labels		
TOS220	Integrity™ V500 Carrying Case		

*XP and G1 Integrity™ V500 Systems have been discontinued. If your VivoLink™ has a serial number beginning with VI or VK, please contact support@vivosonic.com. (Note: G1 and G2 hardware components are not interchangeable.)

Warranty Terms:

All new Integrity™ V500 systems include a ONE-year Standard Equipment Warranty. Warranty Limitations: Disposables, damage due to misuse, system performance issues due to the installation of third party software, re-calibration, equipment shipping & handling cost from installation site to Vivosonic are not covered by warranty. **120-Day warranty limit on battery packs. Laptop computer warranty handled through computer OEM manufacturer.